



## COMBINED MEG/EEG-SOURCE ANALYSIS



- Start Curry: *double-click* the **Curry 7** icon.
- Open the following database:  
**C:\Curry 7 Course\Source Analysis\MEG EEG Somatosensory\MEG EEG.mdb**
- Expand Subject **SEP & SEF**.
- Open study **MEG EGG MRI**.  
(the data import wizard does not appear as data have already been parameterized)
- Select **Butterfly Plot**  (*Alt+B*).
- In **Noise Estimation**, use **Pretrigger Interval**.
- Use **14...26ms** as the analysis timerange (*Ctrl/Alt+left/right cursor keys*).






For combined EEG and MEG analysis, data from both modalities must be stored in the same data file or in separate matching (sampling rate, number of samples) files with MEG first.

## NOTES

## PCA ANALYSIS

Use PCA for MEG and EEG alone and combined to learn about data complexity:

- *Ctrl-click* to change the display but not the parameters.
- From the toolbar, activate a **PCA** .
- In the **Functional Data** parameters, expand **Channel Groups / Rereferencing**.
- Take a look at the PCA components and their SNRs for MEG&EEG combined, MEG alone, and EEG alone:
- Uncheck **EEG**.
- Check **EEG**, uncheck **MEG**.
- Check **EEG** and **MEG**.




MEG data shows higher SNR and less complexity.

## NOTES

## DIPOLE ANALYSIS

Perform dipole fits for the different modalities and observe the achieved goodness-of-fit:

- Switch to the **Maps, 3D View** display.
- In **Dipole Fit**, change the **Dipole Type** to **Moving**.
- From the toolbar, activate **MGFP**, **Dipole Strengths**, and **Goodness-Of-Fit**. The achieved goodness-of-fit is well below the expected goodness-of-fit (dotted line).
- In **Channel Groups / Rereferencing** (click **Functional Data** parameters), uncheck **EEG**. For most samples, achieved and expected goodness-of-fit match.
- Check **EEG** and uncheck **MEG**. Achieved and expected goodness-of-fit match.
- Check both **EEG** and **MEG**.
- In **Head Model** (click **Source Reconstruction** parameters), press **Conductivity Factor: Fit**. The conductivity factor (2.22) is fitted. Achieved and expected goodness-of-fit match.
- Change dipole type **Fixed MUSIC**, set **Number of Dipoles** to **3**.







The MEG-EEG **head model conductivity factor** is a calibration factor for the conductivities which are used for an EEG analysis. (MEG signal strengths only depend on relative conductivities, while EEG signal strengths depend on absolute conductivities.)


*M Fuchs, M Wagner et al., 1998. Improving source reconstructions by combining bioelectric and biomagnetic data. Electroenceph clin Neurophysiol 107:93-111*

## NOTES

## IMAGE ANALYSIS

Load image data and view the dipole results in their anatomical context:

- Switch to the **Image Data** display.
- Make sure **Time Range mode** is activated (press  toolbar button).
- Step through the dipole results using the **Previous/Next Dipole** toolbar buttons  .
- Zoom in and out using the + and – keys.
- Press the **Setup BEM Geometry** toolbar button .
- Change **Create** to **Cortex** and press **Start**.




The **Previous/Next Dipole** toolbar buttons iterate through the active dipole results.

## NOTES

## CUTPLANE OVERLAY

View the dipole results in their anatomical context:

- Switch the display to **3D View** and switch to **Right View** , slightly rotate the cortex surface and zoom in using the mouse wheel.
- Switch off ☐ **Functional Landmarks**, ☐ **Anatomical Landmarks**, ☐ **Electrodes**, and ☐ **Coils**.
- In the ☒ **Cortex (129) 2.5mm** properties set **Cutplane Mode** to **Triple** and **Through** to **3D Cursor**.
- Make sure **Time Range mode** is activated (press toolbar button).
- Step through the dipole results using the **Previous/Next Dipole** toolbar buttons .



- Switch on ☒ **Functional Landmarks**, ☒ **Anatomical Landmarks**, ☒ **Electrodes**, and ☒ **Coils**.
- In the ☒ **Cortex (129) 2.5mm** properties set **Cutplane Mode** to **Off**.





The **Image Data** display can also be used to set the 3D Cursor position for the cutplanes.


## NOTES



## CURRENT DENSITY ANALYSIS

Perform a cortical current density analysis:

- Switch to the  **Maps, 3D View** display, **Top View** .
- Perform an **sLORETA** analysis using the cortical mesh (**Cortex (129) 2.5mm**) as the source space.
- Activate **ICA** and perform a **Component-Based sLORETA CDR** for the first three ICA components.  
(*Ctrl-click* a component to deselect all others)






**Component-Based CDR** allows to compute the source image per ICA component.



## NOTES


## DISPLAY CONTROL

View results within Curry:

- In **Output**, in the upper right corner, press the **Maximize** button . The output field fills the whole display area.
- Press the **Shrink** button .



- Switch the display to **Localize**.
- **Right-click** and select **Import Dipoles**.
- In **Localize**, in the upper left of the location list, press the **Maximize** button . The Localize list is completely visible.
- Press the **Shrink** button .






**Maximize** buttons make it possible to focus on a part of the display or the user interface.

## NOTES

## REPORTS

Paste results and images to the Report:

- Switch to the  **Report**.
- From the **Source Results** menu, select **Report** and **Append Dipole Description**. (or press the  toolbar button). Dipole results are appended to the report.
- In the 3D View display, *right-click* and select **Hardcopy** and **Append Image to Report** (or *Ctrl+Shift+R*). The image is appended to the report.
- *Right-click* the **Report** and select **Open in Editor....** Microsoft Word or another Rich Text Format (rtf) editor opens.
- From the **Edit** menu, select **Options...** and switch to the **Hardcopies** page. On this page, the resolution of hardcopies and movie files can be adjusted.



The **Report** module can be used to collect textual and pictorial results. It is a convenient alternative to saving multiple screenshots.

## NOTES



## RESULTS EXPORT


Export results in Excel and MATLAB formats:

- From the **Source Results** menu, select **Export Results** and **Export Dipoles to Excel...**
- Enter a filename and check **Open in Excel**. Press **Save**.  
A .csv file is written, Excel opens, and dipoles are displayed as a spreadsheet.

A1	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
t [ms]	x [mm]	y [mm]	z [mm]	j [uAmm]	nx	ny	nz	ellipsoid	a [mm]	b [mm]	c [mm]	res dev	nc	var norm	var orig [%]
1	14	-11.7	-29.1	11.4	3.36	0.09	-0.71	-0.7	286.7	26.7	36.2	70.9	44.4	62.7	80.26
2	15	-0.3	-29.1	13.5	4.37	0.07	-0.81	-0.58	68.4	19	22.5	38.1	31.1	92.9	90.35
3	16	8.9	-25.5	18.6	4.4	0.03	-0.8	-0.6	42.4	15	18.9	35.7	31.4	70.9	90.17
4	17	24.8	8.1	61.4	1.41	0.21	-0.91	-0.36	44	16.3	18.5	34.9	43.9	36.7	80.7
5	18	41.7	15.5	84.2	1.3	0.36	-0.93	0.12	1.9	4.4	9.2	10.9	16.9	10	97.15
6	19	41.7	16.4	83.8	3.21	0.3	-0.94	0.18	0.1	1.8	3.7	4.5	14.6	5.32	97.86
7	20	42	16.8	83.6	5.91	0.2	-0.97	0.14	0	1	2	2.3	8.77	4.16	99.23
8	21	42.6	17	83.9	7.23	0.06	-1	0.07	0	0.8	1.6	1.8	7.4	3.49	99.45
9	22	42.8	16.7	83	6.47	-0.22	-0.96	-0.16	0	1	1.9	2.4	10	4.39	98.99
10	23	43.8	17.3	83.6	4.68	-0.53	-0.71	-0.47	0.1	1.9	2.5	4.2	20.1	11.3	95.94
11	24	42.2	17.4	86.3	4.03	-0.63	-0.38	-0.68	0.2	2.7	3.6	5.7	30.6	46.7	90.64
12	25	36.6	20.3	97.5	3.04	-0.55	-0.1	-0.83	0.4	2.8	4.3	8.8	30.4	59.7	90.76
13	26	31.7	16.7	87.1	2.72	-0.56	0.3	-0.77	1	4.5	5.1	10.8	24.2	21.2	94.13
14															95.52

- From the **Source Results** menu, select **Export Results** and **Export Currents to MATLAB...**
- Enter a filename and check **Open in MATLAB**. Press **Save**.  
MATLAB opens (better if already running) and currents are displayed.

Name	Value	Min	Max
V	<18688x4 double>	0	69.4777
curriyedr	<4x18688x13 double>	<Too ...>	<Too ...>
curriydev	<2x13 double>	0.1231	1
curriyloc	<3x18688 double>	-72.0503	106.84...
curriytri	<3x27468 int32>	1	18688
hpatch	174.0013	174.0013	174.00...
nLoc	18688	18688	18688
nflow	4	4	4
nSam	1	1	1
nVal	18688	18688	18688



- Exit Curry



When exporting to **Excel**, a Comma Separated Values (csv) file is written (can be opened in all spreadsheet programs).

When exporting to **MATLAB**, a .mat file is written. If **Open in MATLAB** is checked, a sample .m-file is additionally written, showing how to access the exported data.



## NOTES